


ORDERING INFORMATION

Part Number	Description	Order Qty.
<b>MINILOK</b>		
212850	MINILOK QUICKANCHOR Plus w/ #0 ETHIBOND® (OS-2), 2mm Drill Bit	1/box
212851	MINILOK QUICKANCHOR Plus w/ 2/0 ETHIBOND (V-5), 2mm Drill Bit	1/box
212852	MINILOK QUICKANCHOR Plus w/ 2/0 PANACRYL® (RB-1), 2mm Drill Bit	1/box
<b>MICROFIX</b>		
212855	MICROFIX QUICKANCHOR Plus w/ 3/0 ETHIBOND (V-4), 1.3mm Drill Bit	1/box
212856	MICROFIX QUICKANCHOR Plus w/ 4/0 ETHIBOND (C-1), 1.3mm Drill Bit	1/box
212857	MICROFIX QUICKANCHOR Plus w/ 4/0 ETHIBOND (P-3), 1.3mm Drill Bit	1/box
<b>Mini Anchor</b>		
212033	Mini QUICKANCHOR Plus w/ 2/0 ETHIBOND Excel (V-5)	1/box
212133	Mini QUICKANCHOR Plus w/ #0 ETHIBOND Excel (OS-2)	1/box
211083	Mini Drill, 2.1mm	1/box
213072	Mini Drill Guide	1/box
<b>Micro Anchor</b>		
212865	Micro QUICKANCHOR Plus w/ 3/0 ETHIBOND Excel (V-4), 1.3mm Drill Bit	1/box
212866	Micro QUICKANCHOR Plus w/ 4/0 ETHIBOND Excel (C-1), 1.3mm Drill Bit	1/box
212867	Micro QUICKANCHOR Plus w/ 4/0 ETHIBOND Excel (P-3), 1.3mm Drill Bit	1/box
<b>TACIT Anchor &amp; Instrumentation</b>		
212036	TACIT QUICKANCHOR w/ 1.7mm Hand Drill Bit	1/box
212236	TACIT QUICKANCHOR w/ 1.7mm Cranio Drill Bit	1/box
222220	TACIT Threader (Double End)	2/box
<b>Strada Small Joint Instrumentation</b>		
224200	Grasper, 2.5mm x 90mm	1/box
224201	Straight Punch, Blunt, 2.1mm x 90mm	1/box
224202	Right Punch, Blunt, 2.1mm x 90mm	1/box
224203	Left Punch, Blunt, 2.1mm x 90mm	1/box
224204	Down Punch, Blunt, 2.1mm x 90mm	1/box
224205	Large Punch, 2.7mm x 90mm	1/box
224206	Straight Probe, Marked, 1mm x 65mm	1/box
224207	30° Probe, 1mm x 65mm	1/box
224208	90° Probe, Marked, 1mm x 65mm	1/box
224209	Scissors, 2.5mm x 90mm	1/box
224210	Small Joint Replacement Shear Pin	1/box
215201	Small Joint Arthroscopy Tray	1/box
<b>GII QUICKANCHOR Plus</b>		
222983	GII QUICKANCHOR Plus w/ ORTHOCORD™	1/box
212034	GII QUICKANCHOR Plus w/ ETHIBOND	1/box
212134	GII QUICKANCHOR Plus w/ PANACRYL	1/box
<b>VAPR Family of Electrodes</b>		
225211	2.3mm Short Side Effect Electrode, 85mm	1/box
225213	2.3mm Short Wedge Electrode, 85mm	1/box
225252	2.3mm VAPR TC (Temperature Control) Electrode	1/box


Surgeon's Signature

Date

# Small Bone Solutions




# Small Bone Solutions & Craniofacial Repairs




## MINILOK™ Absorbable QUICKANCHOR™ Plus

MINILOK provides strong, absorbable soft tissue fixation for wrist, hand, foot and ankle applications.




## MICROFIX™ Absorbable QUICKANCHOR Plus

MICROFIX is an excellent solution for precise soft tissue-to-bone fixation in the smaller bones of the digits.




## Mini QUICKANCHOR Plus

This small anchor is for use in certain shoulder, elbow, wrist, hand, foot and ankle applications.




## Micro QUICKANCHOR Plus

This tiny anchor is for use in certain hand and craniofacial applications.




## GII® QUICKANCHOR Plus

The GII QUICKANCHOR Plus is for use in a range of repairs to the shoulder, wrist, hand, foot, ankle and elbow.



## TACIT®

This small anchor is for use in certain wrist, hand and craniofacial applications.



## VAPR® Small Joint Electrodes

The 2.3mm diameter, malleable tip, bipolarity, and short shaft design make these VAPR electrodes useful for hard-to-reach areas of the wrist, ankle and elbow. The solid tip design delivers controlled energy to precisely where it is needed to cut and contour soft tissue.

## Strada Small Joint Instrumentation

Graspers, biters, scissors and probes for small joint arthroscopy applications

