

# Durable Precision


at Your Fingertips

## Conquest Series Manual Instruments

Designed for a broad range of applications, Conquest Instruments combine world-class craftsmanship with unique technical precision. Stryker offers a complete line of 40 instruments including, scissors, graspers, punches and Small Joint Instruments.

### Stryker's Conquest Manual Instruments:

- Ergonomically designed to decrease fatigue.
- Milled from a solid piece of steel and coated with chromium for extra protection.



# Reliable Accuracy


and Control

## Champion Small Joint Instruments

The Champion Small Joint manual instruments are the ultimate tools for small joint arthroscopy. This durable instrument set includes small joint specialty punches, graspers, and scissors that range in size from 2.1mm to 2.7mm.

### Stryker's Champion Small Joint Instruments:

- Set includes up, down, right, and left biters for hard to reach tissue pathology.
- All instruments have a replaceable safety pin mechanism built into the handle. If the instrument experiences an abnormally high load of force, the safety pin will break instead of the working instrument tip.



### Small Joint Probes & Rigid Instruments:

A wide variety of mini specialty tools are available.

