



MEDICAL PRODUCTS

# Shoulder Arthroscopy Instruments

## Compatible for Tape Suture

The next generation of instruments for passing tape suture

### ■ Tape Cutter


An optimized tape suture cutter that cuts 2 strands simultaneously, designed to cut tape up to 2.5mm

### ■ Penetrating Grasper

Shuttling device which passes a suture tape between the jaws through the tissue and releases it in its new position.

### ■ Sixter & Clever Hook

Suture tape passer – allows easy placement and retrieval of suture tape in a rotator cuff procedure.




Instrument designed to fit all types of suture tape from 1.0mm width to 2.5mm width

T.A.G. Medical Products is a leader in the manufacturing of surgical instruments, including high-quality, innovative medical solutions for arthroscopic procedures.


We  Care

# Tape Suture Grasper / Passer


## Mechanical Dimensions




Penetrating  
Grasper 45° Left




Penetrating  
Grasper 45° Right




Penetrating  
Grasper Straight



Penetrating  
Grasper 35° Up



Sixter Left -  
Clever Hook Left



Sixter Right -  
Clever Hook Right

## Tape Cutter

- The Suture Cutter is designed to cut tape up to 2.5mm
- The Suture Cutter is designed to cut 2 strands of tape suture simultaneously
- The Suture Cutter is designed for cutting without visualization
- The Suture Cutter includes a ratchet mechanism to prevent premature tape cutting
- The Suture Cutter has capability of knot pushing & preventing knot cutting

**tag** arthro



**T.A.G. Medical Products Corporation Ltd**

Kibbutz Gaaton 2513000 Israel | Tel: +972-(0)4-9858400 | [www.tag-med.com](http://www.tag-med.com)